



**Centro Tessile Serico Spa consortile**

Via Castelnuovo, 3 - 22100 Como, Italia

Tel +39 031 331211/120 Fax +39 031 3312180

Amministrazione/Contabilità +39 031 3312102-154

P.IVA - C.F. 02381270137

Capitale sociale euro 552.693,18 i.v.

mailbox@textilecomo.com

www.textilecomo.com

## Material quantity for tests

Below are the general indications regarding the minimum quantities of material necessary for the execution of the main tests. Except for agreements between the parties for special cases, the quantity indicated guarantees the execution of the tests in a complete and correct manner.

TEST ON YARNS	1 cone
DEFECTS AND EVENNESS OF RAW SILK	12 cones with at least 100 grams of yarns
PHYSICAL STRUCTURAL TESTS ON TEXTILES	1 meter total width
CHEMICAL TESTS ON TEXTILES	50 grams
DIMENSIONAL STABILITY TESTS	1 meter total width
COLOURFASTNESS TESTS	1 cone or 20 cm total width (all colours inclusive)
CHEMICAL PRODUCTS AND AUXILIARY ANALYSIS	100 grams
WASTEWATER ANALYSIS	2 liters
FLAMMABILITY TESTS	
Curtains - Italy	6 m <sup>2</sup>
Curtains - Foreign	3 m <sup>2</sup>
Upholstery / Mattresses:	
- Covering	3 m <sup>2</sup>
- Padding	to be agreed
Apparel – USA	1 m <sup>2</sup>
Protective apparel	2 m <sup>2</sup>
Interior vehicles	1 m <sup>2</sup>



**Centro Tessile Serico Spa consortile**

Via Castelnuovo, 3 - 22100 Como, Italia

Tel +39 031 331211/120 Fax +39 031 3312180

Amministrazione/Contabilità +39 031 3312102-154

P.IVA - C.F. 02381270137

Capitale sociale euro 552.693,18 i.v.

mailbox@textilecomo.com

www.textilecomo.com

## Samples storage methods

The following are the indications to be applied by the customer before sending or delivery the sample to the laboratory

### FREE AND HYDROLYZED FORMALDEHYDE

*Deliver the samples separately in a closed polyethylene bag, wrapped in aluminium foil*

### MOISTURE CONTENT AND/OR MOISTURE REGAIN

*Deliver well sealed sample*

### TESTS OF ABRASION / PILLING

*Deliver minimum patten repeat*

### SURFACE REPLICA

*Mark the defect*

### DETERMINATION TO WATER PENETRATION - CONSTANT PRESSURE TEST AND CREASY RECOVERY TEST

*Important not crease the sample*

### ELECTROSTATIC PROPERTIES – SURFACE RESISTANCE AND ELECTRICAL RESISTANCE

*Deliver in a closed paper envelope*

### WASTEWATER ANALYSIS

*For collection and transport of water samples use clean glass or polyethylene containers, which must be kept refrigerated ( $4 \div 10$  ° C).*

The Laboratory Secretary and the Customer Service are available  
for information and clarifications

[segreteria@textilecomo.com](mailto:segreteria@textilecomo.com)

+39 031 3312121 - 170

[scienti@textilecomo.com](mailto:scienti@textilecomo.com)

+39 031 3312122